

To be completed at Charles River



charles river

LTM Order Number:

Receipt Number:

Date Received:

Cell Line and Research Biologics PCR Testing Submission Form

For PCR Rodent Infectious Agent (PRIA), Direct Sampling, and Exhaust Air Dust (EAD) submissions, please use the [PRIA PCR Submission Form](#).

LTM Customer ID#* _____

New customers, write NEW in this field. Existing customers can find their LTM customer ID# in the header of result reports from previous submissions. Need help? Contact LabServices@crl.com.

Customer Information

* = required

Bill To Information

Same as Customer

Payment information required at the time of submission.

Institution* _____

Institution* _____

Pay by PO

Address 1* _____

Address 1* _____

PO # _____

Address 2 _____

Address 2 _____

Pay by Credit Card (CC)

City* _____

City* _____

Contact _____

State* _____ Zip Code* _____

State* _____ Zip Code* _____

Phone _____

Contact* _____

Contact* _____

To pay with a CC on file, indicate the last four (4) digits and expiration date. For CC security, do not write the full CC number.

Phone* _____

Phone* _____

CC Last 4 Digits _____

Email* _____

Email* _____

CC Exp Date _____

Sample Information

To the best of my knowledge these samples and/or specimens do not contain any infectious agent or material which might pose a threat to human health. If this is not the case, please call Charles River before sending samples.

Signed: _____

Date: _____

Total Number of Samples: _____

Ship Date: _____

Duplicate samples included for retesting Custom panel(s) requested; panel name(s) entered below Add Contamination CLEAR Panel

Vial #	Sample ID	Service Requested	Sample Type	Species
1				
2				
3				
4				
5				
6				
7				
8				
9				
10				

Notes: _____

Instructions to lab: _____

Want an easier way to fill out a submission form and track your testing with Charles River? Try using LTM™ -- the online order wizard will guide you through your testing request. Visit www.criver.com/lrm to get started.



Shipping Instructions

Please only print and include page 1 with your samples.

Questions? Email LabServices@crl.com or visit www.criver.com/info/dx

Infectious Disease PCR Testing is available for research purposes only.

Sample Specifications:

Tissues: Place entire tissue sample in a 2.0 mL microcentrifuge tube. Tissue sample may represent one animal, or poled animals, as long as the total size of tissue is not more than approximately a 500 μ L volume.

Serum or Urine: Send a maximum of 500 μ L undiluted sample.

Research Biologics:

Cell lines: Send two (2) undiluted 200 μ L samples of each line; the second vial allows for confirmation of positives. Cell number is not critical; however, please note in the "Instructions to lab" section on page 1 if there are more than 5×10^7 cells/mL. If possible, samples for Mycoplasma testing should be passaged without antibiotics prior to shipping. Cells do not need to be viable for PCR testing. If shipping samples containing DMSO, please indicate this in the "Instructions to lab" section on page 1 to facilitate proper disposal.

Antibodies and other biologics: Send two (2) undiluted 200 μ L samples. If these volumes are not possible, please contact us for recommendations. Please indicate if the protein concentration is greater than 1.5 mg/mL. If possible, samples for Mycoplasma testing should be passaged without antibiotics prior to shipping.

Storage and Sample Submission:

Store frozen; ship overnight on dry ice.

Place a copy of page 1 of this form in the appropriate shipping container with your samples. Please ship following International Air Transport Association (IATA) guidelines; visit www.criver.com/shippingmaterials to request free sample collection and shipping materials online.

Submit your samples to Charles River:

**Charles River
RADS
PCR Laboratory
261 Ballardvale Street
Building 22 Loading Dock
Wilmington, MA 01887
USA**

**Charles River
RADS
West Coast Laboratory
10788 Roselle Street
San Diego, CA 92121
USA**

**Charles River
RADS
327 impasse du Domaine Rozier
Domaine des Oncins
69210 SAINT GERMAIN NUELLES
FRANCE**